

Powerheart G3 Elite Storage & Signage Options



Carry Bag Solutions

G3 Premium Carry Case

Part No. XCAAED001A

This semi-rigid carry bag provides additional protection for your Powerheart G3 AED and makes it easy to grab the unit in a hurry. The large back pouch is perfect for storing extra electrodes and the shoulder strap is easily adjustable. In addition, the Rescue Ready indicator and pads' expiration date are visible at all times. Fits all 9300 series AEDs.



Pelican Case

Part No. XCAAED003A

This hard-sided, waterproof case secures the Powerheart G3 in rugged environments, and makes it easy to transport with reliable protection from bumps, drops, and the elements.



Signage and Storage

3D Wall Sign Kit

Part No. 168-0056

This 3D wall sign kit, helps ensure First Responders know where to locate your AED in a sudden cardiac arrest emergency. The kit includes a highly visible sign which projects from the wall and a decal sticker.



AED Wall Hanger / Bracket

Part No. 180-2023

The AED Wall Hanger / Bracket together with signage is an economical way to store your AED on the wall and within easy access. The wall hanger is suitable for the AED by itself or if the unit is inside a carry case.



AED Wall Storage Sleeve

Part No. XWSAED001A

This storage sleeve is an economical, attractive way to keep your AED in an easily accessible location and ensure high visibility to rescuers. The storage sleeve with security belt is tested to the EN 1789 standard for use in emergency vehicles. This Metal Sleeve is for use with the G3 Premium Carry Case.



Wall Cabinets

AED Wall Cabinet with Alarm

Part No. 55-00392-30

The distinctive curved shape and green door makes it easy to find and open with AED signage. Wall-mountable steel AED (Automated External Defibrillator) cabinet is specifically designed to securely store your defibrillator in its case for easy access. The cabinet comes with a battery alarm which sounds when the door is opened. The cabinet has the option to be locked with a key or kept unlocked. The alarm can easily be disabled or enabled.



AED Storage Cabinet with Strobe Light and Alarm

Part No. 50-00392-30

These well-constructed AED wall cabinets store your AED in an easily accessible location. The AED sits inside the white cabinet with a see-through door.

Security options: Make it widely-known that an emergency is in progress. These options are triggered when the door is opened. The alarm can be disabled and enabled using a hidden internal switch.



AED Waterproof Cabinet

Part No. XWCAED001

Our waterproof cabinet ensures your AED is safely stored and protected from any water damage or harmful liquids. The cabinet is suitable for all environments, particularly aquatic centres and surf-living saving clubs.

