



FRED® easy® and FRED® easyport® Trainer

Train yourself to save lives!

Currently more and more responders of the chain of survival equip themselves with automated external defibrillator: First aid associations, companies, swimming pools, public places and cities introduce general public defibrillation programs.

All these initiatives express a more and more important training need, that's why SCHILLER developed its range of training devices:

FRED® easy® Trainer and FRED® easyport® Trainer are training devices designed to teach users of AED's how to save lives in a realistic way.

Features:

- Operation identical to the real devices, display with instructions and vocal messages
- 9 preset scenarios
- Infrared remote control, enabling the instructor to choose and control a scenario
- Training electrodes, suitable for any mannequin
- Reusable cable



FRED® easy® Trainer

FRED® easy® Trainer is an easy-to-handle training device to teach AED users how to save lives.

2 in 1 defibrillator: fully automatic and half automatic to be able to train on all types of devices.

Specifications of FRED® easy® Trainer:

- Configurable in semi-automatic or fully automatic mode
- Working on cells or rechargeable battery



FRED® easy® Trainer

FRED® easyport® Trainer

FRED® easyport® Trainer is the world's smallest training device to teach AED users how to save lives.

Specifications of FRED® easyport® Trainer:

- Semi-automatic mode
- Power supply: 9 V battery



FRED® easyport® Trainer



Technical data:

FRED® easy® Trainer

Dimensions:

230 x 220 x 70 mm (l x w x h)

Weight:

< 1 kg

Power supply:

- Alkaline batteries (4x1.5 V)
- Ni-Cd battery (6 V): to be ordered separately (P/N 4-07-0010)
- Battery charger: to be ordered separately (P/N 3-55-0030)

Accessories:

Delivered with one remote control, one cable for training electrodes, 10 sets of training electrodes, alkaline batteries (4x1.5 V).

Technical data and pictures are subject to change without notice.

Technical data:

FRED® easyport® Trainer

Dimensions:

133 x 126 x 35 mm (l x w x h)

Weight:

< 300 g

Power supply:

Alkaline battery (9 V)

Accessories:

Delivered with one remote control, one cable for training electrodes, 10 sets of training electrodes, alkaline battery (9 V)

Technical data and pictures are subject to change without notice.

Asia

SCHILLER Asia-Pacific

7 Jalan Ara SD 7/3B, Bandar Sri Damansara
52200 Kuala Lumpur, Malaysia
Phone +603 6272 3033
Fax +603 6272 2030
sales@schiller.com.my
www.schiller-asia.com

Austria

SCHILLER Handelsgesellschaft mbH

Kampfmüllerweg 24
A-4040 Linz
Phone +43 732 70 99-0
Fax +43 732 757 000
sales@schiller.at
www.schiller-austria.com

C.I.S. & Russia

SCHILLER AG

Moscow Representative Office
Office 403, 1-st Yamskogo Polyta 15
125124 Moscow, Russia
Phone +7 (495) 970 11 33
Fax +7 (495) 956 29 10
mail@schiller-ag.com
www.schiller-cis.com

France

SCHILLER Médical S.A.S.

4, rue Louis Pasteur
F-67162 Wissembourg/Cedex
Phone +33 3 88 63 36 00
Fax +33 3 88 94 12 82
info@schiller.fr
www.schiller-france.com

France (office Paris)

SCHILLER Médical S.A.S.

7, rue Raul Follereau, Bussy Saint Georges
F-77608 Marne la Vallée/Cedex
Phone +33 1 64 66 50 00
Fax +33 1 66 50 10
info@schiller.fr
www.schiller-france.com

Germany

SCHILLER Medizintechnik GmbH

Rudolf-Diesel-Strasse 14
D-85521 Ottobrunn
Phone +49 89 62 99 81-0
Fax +49 89 609 50 90
info@schillermed.de
www.schiller-germany.com

India

SCHILLER Healthcare India Pvt. Ltd.

D.C. Silk Mills Compound, 'A' Wing, 1st floor
5, Chunawala Estate, Kondivita Lane
Andheri-Kurla Road, Andheri (E)
Mumbai - 400 059, India
Phone +91 22 2826 35 20
Fax +91 22 2826 35 25
sales@schillerindia.com
www.schillerindia.com

Japan

SCHILLER Japan, Ltd.

Biomedical Research Building 218
1-2-3 Kasumi, Minami-ku
Hiroshima 734-8551
Phone +81 82 250 2055
Fax +81 82 253 1713
koji.maekawa@schiller.jp
www.schiller.jp

Spain

SCHILLER ESPAÑA, S.A.

Ctra. de la Coruña
A6 - km 18,200
E-28230-Las Rozas/Madrid
Phone +34 91 713 01 76
Fax +34 91 355 79 33
schiller@schiller.es
www.schiller-spain.com

Switzerland

SCHILLER-Reomed AG

Riedstrasse 14
CH-8953 Dietikon
Phone +41 44 744 30 00
Fax +41 44 740 37 10
sales@schiller-reomed.ch
www.schiller-reomed.com

Turkey

SCHILLER TÜRKİYE

Haliç Rifat Paşa Mah., Darülaceze Cadd.
Perpa Ticaret Merkezi, B/Blok No: 41 K: 2
Okmeydanı-Sisli - İstanbul
Phone +90 212 210 8681 (pbx)
Fax +90 212 210 8684
sales@schiller-turkiye.com
www.schiller-turkiye.com

USA

SCHILLER America Inc.

2131 N.W. 79th Avenue
Doral, Florida 33122
Phone +1 786 845 0620
Fax +1 786 845 0602
sales@schilleramerica.com
www.schilleramerica.com

08/07




No d'article : 2.500496



CARDIAC
defibrillators

www.cardiacdefibrillators.com.au
sales@cardiacdefibrillators.com.au
(03) 9429 2666



SCHILLER
The Art of Diagnostics